



Amendments to the Specification

Please amend the "Cross Reference to Related Applications" section, appearing on page 1 of the specification, as follows:

This Application is a continuation of U.S. Application No. 09/984,664, filed October 30, 2001, ~~and published on May 29, 2003 as patent application publication number US 2003-0099950~~ now U.S. Patent No. 7,045,319 B2, issued May 16, 2006; which ~~are~~ is incorporated by reference herein in ~~their~~ its entirety.

Please amend paragraph [0051] as follows:

[0051] FIG. 13: Deamination conversion of unmethylated cytosine groups in DNA. Deamination converts unmethylated C to U. Methylated C groups, such as those in CpG islands that regulate eukaryotic genes, are resistant to deamination and remain as C in the product DNA. If 100% deamination occurs, methylated DNA will still contain CpG doublets, whereas unmethylated DNA will contain no cytosine and will now contain UpG where CpG doublets were before deamination. This difference in DNA sequence can be used to distinguish between methylated and unmethylated DNA by abortive transcription because the two DNAs encode different dinucleotides. (SEQ ID NOS: 1, 2, 3 and 17).

Please amend paragraph [0052] as follows:

[0052] FIG. 14: Detection of methylation using dinucleotide synthesis. Dinucleotide synthesis can be used to assess the overall methylation state of DNA. In the presence of RNA polymerase, CTP or a CTP analog (R1-C-OH), and GTP or a GTP analog (R1-CpG-R2), the deaminated methylated DNA template will produce n copies of a labeled dinucleotide product, where n is proportional to the number of methylated CpG

dinucleotides in the starting DNA. The deaminated unmethylated DNA template can produce no dinucleotide with these substrates because the template no longer encodes "C" at any position. (SEQ ID NOS: 4 and 5).

Please amend paragraph [0055] as follows:

[0055] FIG. 15: Assessing methylation status of specific CpG sites in CpG islands by abortive initiation. Target site probes can be used to examine the methylation status of specific CpG islands in specific genes. In the deaminated methylated DNA, the dinucleotide CpG is encoded by the template at the 3 methylated sites 1, 3 and 4, but not by the unmethylated site 2. To specifically determine if Site 3 is methylated and if so, to what extent, position (C21) can be targeted with a Target Site Probe, as described in Figure 11. The template C in question is positioned at the junction of the bubble region and the downstream duplex so that it encodes the next incorporated nucleotide for appropriately primed RNA polymerase that binds to the bubble region. If a labeled initiator R_1-N_XpC-OH is used, where N_X may be C for a dinucleotide CpC initiator or N_X may be CpC for a trinucleotide initiator, etc., the initiator can be elongated with a labeled GTP analog $pppG-R_{2G}$ to form a trinucleotide $R_1N_XCpG-R_2$. Similarly, if the C in question was not methylated, the position will now be a U and will encode nucleotide A. If an ATP analog $pppA-R_{2A}$ is present, it will be incorporated opposite positions where the C was not methylated. If the GTP analog is labeled with group R_{2G} , which is an energy acceptor from the R group on the initiator, R_1 , then the amount of $R_1N_XCpGR_{2G}$, which will be proportional to the amount of methylated C present at that position, can be quantified by measuring the emission from R_{2G} at wavelength λ_{2GE} . The similar situation exists for incorporation of the ATP analog and measurement of the emission from its R

group, also an energy acceptor from the initiator R_1 . By determining the ratio of the magnitude of emission from the GTP analog to the total emission from both the ATP and GTP analogs, the site can be assigned a methylation index M . If all of the Cs at that position are methylated, $M = 1$. If none of the site is methylated, $M = 0$. (SEQ ID NOS: 6, 7 and 8).

Please amend paragraph [0059] as follows:

[0059] FIG. 19: Signal Generation from abortive promoter. An Abortive Promoter Cassette (APC) consists of one or more oligonucleotides or polynucleotides that together create a specific binding site for an RNA polymerase coupled to a linker region (APC linker) for attachment to target molecules (DNA, RNA, Protein). The APC may contain an artificial promoter, or it may contain the promoter for a specific RNA polymerase. For example, trinucleotide or tetranucleotide products that could be generated from with a common phage RNA polymerase can be made with a labeled GpA or GpApA initiator and a labeled pppG or pppA terminator. (SEQ ID NOS: 9 and 10).

Please amend paragraph [0065] as follows:

[0065] FIG 25. Detection of telomerase activity with reiterative oligonucleotide synthesis. Reiterative oligonucleotide synthesis with DNA polymerases can also be used for signal generation, however, the product oligonucleotides need not be released, but may be joined tandemly in the product. As an example, telomerase activity can be detected by immobilizing a telomerase-specific probe to a solid matrix to capture cellular telomerase, which carries its own RNA template for DNA synthesis. For example, with human telomerase, the RNA template on the enzyme encodes the DNA sequence GGGTTA. The capture probe may contain the sequence GGGTTA, which will be added

reiteratively to the end of the telomerase capture probe, if telomerase is present in the sample. Signal generation can be achieved in several ways, one of which involves including one or more reporter tagged dNTPs in the synthesis reaction to produce a product that has multiple R_1 groups attached along the backbone of the DNA product. For detection, this product can then be hybridized to a complementary probe containing nucleotides with a second R group (R_2) attached that will hybridize to the R_1 labeled product. This will bring the R_1 and R_2 groups together for signal generation via FRET from between R_1 and R_2 , or via other methods. Alternatively, telomerase may incorporate 2 labeled nucleotides in the product DNA and look for energy transfer between the 2 labeled nucleotides in the single strand of DNA. (SEQ ID NOS: 11, 12, 13 and 14).

Please amend paragraph [0070] as follows:

[0070] FIG. 29A and B. Portion of the contig sequence of the CDKN2A gene. The sequence represents a small portion of the contig starting at 856630 nucleotides from the start of the contig sequence. The sequence represents a CpG island. Contig number: NT_008410.4. (SEQ ID NO: 15).

Please amend paragraph [0247] as follows:

[0247] ATATACTGGGTCTACAAGGTTTAAGTCAACCAGGGATTGAAATATAACTTTT
AAACAGAGCTGG. (SEQ ID NO:16). The DNA sample is incubated with the
capture probe to allow hybridization. A representative hybridization protocol is as
follows: (1) prehybridize with 2.5X SSC, 5X ~~Denhardt's~~ Denhardt's at room temperature
for 30 minutes; (2) hybridize with 2.5X SSC, 5X ~~Denhardt's~~, Denhardt's 30%

formamide at room temperature for 2 hours; (3) wash twice with 1X SSC at 42°C for 10 minutes, maintaining 42°C; and (4) wash three times with 0.1X SSC at 42°C for 10 minutes, maintaining 42°C.